

BIOLOGICAL BRILLOUIN SCATTERING WITH THE HYPERFINE SPECTROMETER

LightMachinery Inc., 80 Colonnade Road North, Ottawa ON Canada
For questions, contact hyperfine@lightmachinery.com

LightMachinery
Excellence in Lasers and Optics


INTRODUCTION

A rapidly growing body of work shows that mechanical cues play a key role in regulating cellular behaviors and tissue homeostasis. Deregulation of such mechanical properties are thought to be essential in the onset and progression of many diseases, including cancer, osteoporosis, and atherosclerosis¹⁻⁴. Characterizing these physical changes constitute a powerful tool to better understand their underlying mechanism, explore novel therapies (e.g. synthetic matrices⁵), and improve early diagnosis (e.g. ocular and dental fields).


LightMachinery has developed the Brillouin Hyperfine spectrometer to optically probe the mechanical properties of biological materials with unparalleled sensitivity. Its high contrast and ultrafast acquisition are indispensable for successful 3D Brillouin microspectroscopy.

HYPERFINE SPECTROMETER OPERATING PRINCIPLE

Schematic layout of the optical components⁷. The etalon disperses the incident light in the vertical direction, providing a high resolution, but also overlapping orders. The grating separates the overlapping orders in the horizontal direction. The Pump Killer (not shown) is a tunable dual pass air-spaced etalon in reflection.





BRILLOUIN HYPERFINE SYSTEM WITH PUMP KILLER




- Brillouin spectrum of methanol, 250 msec, 25 mW at 532 nm.
- Signal-to-noise ratio ~ 200.
- The pump signal is completely suppressed (with 2 Pump Killer modules, > 70 dB suppression is achieved).

- > 120 dB of effective contrast achievable with the full system.
- Ultra-fast acquisition⁶, massive dynamic range.
- Customizable excitation wavelength from UV to IR.
- Customizable resolution from 0.25 GHz to 5 GHz.
- Up to several THz of spectral range captured simultaneously.




BRILLOUIN SCATTERING SPECTRA – PRELIMINARY RESULTS




a. Brillouin spectrum of an apple slice, 2 sec, 5 mW at 660 nm. b. Raw sensor image corresponding to the spectrum.




a. Brillouin spectra of blood (green) and water (blue). Brillouin frequency shift between water and blood ~ 0.4 GHz. 15 mW at 532 nm, blood spectrum: 1 sec (avg x10), water spectrum: 200 msec (avg x10). b. Excitation configuration showing blood wells.



a. Bovine tooth sample on the Brillouin confocal microscope. b. Excitation configuration showing enamel and dentin samples. Samples provided by MSTATT LLC.



- Brillouin spectra of bovine enamel (green) and dentin (red).
- 10 sec (avg x20), 15 mW at 532 nm.
- In collaboration with MSTATT LLC.



- Brillouin spectra of bovine enamel before (blue) and after (green) soaking for 2 hours in Coca Cola. Brillouin shift of enamel ~ 1 GHz.
- 10 sec (avg x20), 15 mW at 532 nm.
- In collaboration with MSTATT LLC.

CONCLUSIONS



LightMachinery is proudly joining this brilliant collective effort to push the boundaries of biological Brillouin scattering technologies. We are driven by the exploration of uncharted waters in fundamental and applied biological physics.

References

1. Robert Prevedel, Alba Diz-Muñoz, Giancarlo Ruocco, and Giuseppe Antonacci. "Brillouin microscopy: an emerging tool for mechanobiology." *Nature methods* 16, no. 10 (2019): 969-977.
2. Karin A. Jansen, Dominique M. Donato, Hayri E. Balciglu, Thomas Schmidt, Erik H.J. Danen, and Gijsje H. Koenderink. "A guide to mechanobiology: where biology and physics meet." *Biochimica et Biophysica Acta (BBA)-Molecular Cell Research* 1853, no. 11 (2015): 3043-3052.
3. Jay D. Humphrey, Eric R. Dufresne, and Martin A. Schwartz. "Mechanotransduction and extracellular matrix homeostasis." *Nature reviews Molecular cell biology* 15, no. 12 (2014): 802-812.
4. Diana E. Jaalouk, and Jan Lammerding. "Mechanotransduction gone awry." *Nature reviews Molecular cell biology* 10, no. 1 (2009): 63-73.
5. Kyle H. Vining, and David J. Mooney. "Mechanical forces direct stem cell behaviour in development and regeneration." *Nature reviews Molecular cell biology* 18, no. 12 (2017): 728-742.
6. Guqi Yan, Antony Bazir, Jérémie Margueritat, and Thomas Dehoux. "Evaluation of commercial Virtually Imaged Phase Array and Fabry-Pérot based Brillouin spectrometers for applications to biology." *Biomedical optics express* (2020).
7. John Reid, Hubert Jean-Ruel, and Greg Dodd. "Introduction to the operating principles of the HyperFine spectrometer." White paper available on LightMachinery's website. <https://lightmachinery.com/media/1857/hyperfine-principles-of-operation.pdf>

BioBrillouin2020
Virtual Meeting
9th – 11th September