





REVISION HISTORY			
REV	DESCRIPTION	DATE	APP'D
A	INITIAL RELEASE		



DIMENSIONS IN [] ARE MILLIMETERS.

UNLESS OTHERWISE SPECIFIED
DIMENSIONS ARE IN INCHES
TOLERANCES ARE:
DECIMALS ANGLES
.X ± .020 .X ± .5°
.XX ± .010 .XX ± .01°
.XXX ± .005



LightMachinery
Excellence in lasers and optics
80 COLONNADE RD.
NEPEAN, ONTARIO, K2E 7L2, CANADA

MATERIAL
AL-6061-T6

FINISH
BLACK ANODIZED

DO NOT SCALE DRAWING

APPROVALS	DATE
DRAWN BEN SHIM CHECKED H. JEAN-RUEL ENG APPR.	23-Jan-2020
MFG APPR.	

TITLE:
**HORNET 700 - 1070nm
25-40pm RESOLUTION**

SIZE DWG. NO. REV
B HN-9353 A

SCALE: 1:2 SHEET 1 OF 1

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS
DRAWING IS THE SOLE PROPERTY OF
LIGHTMACHINERY. ANY REPRODUCTION IN
PART OR AS A WHOLE WITHOUT THE
WRITTEN PERMISSION OF LIGHTMACHINERY
IS PROHIBITED.